

3 DOF HOVER

The 3 DOF Hover system is ideally suited to study control concepts and theories relevant to real world applications of flight dynamics and control in vertical lift off vehicles.

STUDY FLIGHT DYNAMICS AND CONTROL OF VERTICAL LIFT-OFF VEHICLES



The 3 DOF Hover experiment provides an economical test bed to understand and develop control laws for flight dynamics and control of vehicles with vertical lift off.



System specifications on reverse page.

3 DOF HOVER WORKSTATION COMPONENTS

3 DOF Hover plant Q8-USB data acquisition device VoltPAQ-X4 four channel linear voltage amplifier QUARC real-time control software for MATLAB®/Simulink® Laboratory Guide and User Manual (provided in digital format) Sample pre-built controllers and complete dynamic model



3 DOF Hover workstation

HOW IT WORKS

The 3 DOF Hover consists of a planar round frame with four propellers. The frame is mounted on a three degrees of freedom pivot joint that enables the body to rotate about the roll, pitch and yaw axes. The propellers are driven by four DC motors that are mounted at the vertices of the frame.

The propellers generate a lift force that can be used to directly control the pitch and roll angles. The total torque generated by the propeller motors causes the body to move about the yaw axis. Two of the propellers are counter-rotating, so that the total torque in the system is balanced when the thrust of the four propellers is approximately equal.

The voltage signals going to the motors, as well as the pitch and yaw encoder signals are transmitted through a slip ring. The slip ring removes the need for wires and allows for 360 degrees free motion about the yaw axis. Furthermore, it reduces the amount of friction and loading about the moving axis.




Figure 1

An example of a free body diagram for the 3 DOF Hover provided in the courseware. Students can use this model to study forces interacting with the hover and the direction conventions.

made for science

SYSTEM **SPECIFICATIONS 3 DOF Hover**

CURRICULUM TOPICS PROVIDED

• Derivation of simple dynamic model

• State space representation

pitch and yaw axes

DEVICE SPECIFICATIONS

Device height (ground to top of base)

Encoder resolution (in quadrature)

Motor / propeller force-thrust constant

Motor / propeller torque thrust constant

measurements

Device mass

Helicopter body mass

Helicopter body length

Pitch angle range

Propeller diameter

Motor armature resistance

Motor current-torque constant

Propeller pitch

Yaw angle range

Base dimensions – $W \times L$

FEATURES

Propeller shield Motors actuate pitch Encoders used and roll direction to meaure pitch and roll direction Slip ring allows 360 degree Encoder used rotation about to measure yaw axis vaw angle State feedback control • LQR control design Control parameter tuning • Three degrees of freedom (3 DOF) - body rotates about • Easy-connect cable and connectors · Precise, stiff and heavy-duty machined components • Propellers driven by high-quality Pittman DC motors • Fully compatible with MATLAB[®]/Simulink[®] and LabVIEW[™] • High-resolution optical encoders for precise position • Fully documented system model and parameters provided for MATLAB[®]/Simulink[®], LabVIEW[™] and Maple[™] • Slip ring allows infinite motion about the yaw axis • Open architecture design, allows users to design their own controller 3.46 kg 45 cm 1.39 kg 48 cm 17.5cm × 17.5 cm 8192 counts/rev 75 (± 37.5 deg) 360 deg 0.119 N/V 0.0036 N.m/V 20.3 cm 15.2 cm 0.83 Ω

COMPLETE WORKSTATION COMPONENTS

Plant	3 DOF Hover
Control design environment	Quanser QUARC [®] add-on for MATLAB [®] /Simulink [®]
	Quanser Rapid Control Prototyping (RCP) Toolkit add-on for LabVIEW™
Documentation	User Manual and Laboratory Guide
Real-time targets	Microsoft Windows® and NI CompactRIO
Data acquisition devices	Quanser Q8-USB, QPIDe, or equivalent NI DAQ device supported by QUARC
	NI CompactRIO with four Quanser Q1-cRIO modules
Amplifier	Quanser VoltPAQ-X4 linear voltage amplifier

0.0182 N.m/A

About Quanser:

Quanser is the world leader in education and research for real-time control design and implementation. We specialize in outfitting engineering control laboratories to help universities captivate the brightest minds, motivate them to success and produce graduates with industry-relevant skills. Universities worldwide implement Quanser's open architecture control solutions, industry-relevant curriculum and cutting-edge work stations to teach Introductory, Intermediate or Advanced controls to students in Electrical, Mechanical, Mechatronics, Robotics, Aerospace, Civil, and various other engineering disciplines.

madelforscience

made-for-science.com | quanser@made-for-science.com | +49 69 2474906 10

Products and/or services pictured and referred to herein and their accompanying specifications may be subject to change without notice. Products and/or services mentioned herein are trademarks or registered trademarks of Quanser inc. and/or its affiliates. MATLAB® and Simulink® are registered trademarks of the MathWorks, Inc. LabVIEW[™] is a trademark of National Instruments. Maple[™] is a trademark of Maplesoft. Windows[®] is a registered trademark of Microsoft Corp. ©2015 Quanser Inc. All rights reserved.