

HEAT FLOW EXPERIMENT

The Heat Flow Experiment is a process control plant ideal to demonstrate simple first order model control concepts as well as fluid dynamics and thermodynamics control fundamentals.

A SYSTEM FOR STUDYING SIMPLE DYNAMICS


The Quanser Heat Flow Experiment demonstrates control topics related to fluid dynamics and thermodynamics. It introduces a “bump test” method, a simple technique used to find the first-order model of the system and helps students better understand temperature control strategies, such as an on-off control scheme using a relay switch, and advanced PID control topics including set-point weight and integral anti-windup.

HOW IT WORKS

The Quanser Heat Flow Experiment plant consist of a fiberglass chamber that is equipped with a coil-based heater and a blower at one end and three temperature sensors located equidistantly along the chamber.

The Heat Flow Experiment has a built-in amplifier to deliver power to the heater and blower. The power delivered to the heater and blower is controlled using analog signals. The fan speed can also be controlled using an analog signal, and is measured by a tachometer mounted on the blower measuring angular rate of the blower fan.


The three temperature sensors in the chamber are fast settling platinum transducers.



System specifications on reverse page.

HEAT FLOW EXPERIMENT WORKSTATION COMPONENTS


- Heat Flow Experiment plant
- QB-USB data acquisition device
- QUARC real-time control software for MATLAB®/Simulink®
- Instructor and Student Workbooks, and User Manual (provided in digital format)
- Sample pre-built controllers



The block diagram helps students understand how to implement on-off control using a relay switch.

SYSTEM SPECIFICATIONS

Heat Flow Experiment



CURRICULUM TOPICS PROVIDED

- Experimental modeling – step response
- Transfer function representation
- Model validation
- PID control
- Set-point weighting
- Integrator windup
- Feed-forward

FEATURES

- Fast settling platinum temperature transducers to measure the chamber air temperature
- Tachometer for measuring angular rate of the blower fan
- Coil-based heater
- Solid Plexiglas chamber
- Fully compatible with MATLAB®/Simulink®
- Fully documented system models and parameters provided for MATLAB®, Simulink® and Maple™
- Open architecture design, allowing users to design their own controller

DEVICE SPECIFICATIONS

Device dimensions (L x H x W)	50 cm x 15 cm x 10 cm
Device mass	0.5 kg
Cross sectional chamber area	0.0064 m ²
Blower nominal input voltage	6 V
Blower nominal airflow	1.02 m ³ /min
Maximum wind speed	159.4 m/min
Blower maximum speed	2700 RPM
Tachometer calibration gain	1064 RPM/V
Heater maximum power (at 5 V)	400 W
Temperature sensor calibration gain	20 °C/V
Temperature sensor settling time	4 s

COMPLETE WORKSTATION COMPONENTS

Plant	Heat Flow Experiment
Control design environment	Quanser QUARC® add-on for MATLAB®/Simulink®
Documentation	User Manual, Instructor and Student Workbooks
Real-time targets	Microsoft Windows®
Data acquisition devices	Quanser Q8-USB, QPIdE
The system model and a sample controller(s) are supplied	

About Quanser:

Quanser is the world leader in education and research for real-time control design and implementation. We specialize in outfitting engineering control laboratories to help universities captivate the brightest minds, motivate them to success and produce graduates with industry-relevant skills. Universities worldwide implement Quanser's open architecture control solutions, industry-relevant curriculum and cutting-edge work stations to teach Introductory, Intermediate or Advanced controls to students in Electrical, Mechanical, Mechatronics, Robotics, Aerospace, Civil, and various other engineering disciplines.