



STUDENT WORKBOOK

Ball and Beam Experiment for LabVIEW™ Users

Standardized for ABET* Evaluation Criteria

Developed by:
Jacob Apkarian, Ph.D., Quanser
Paul Karam, B.A.S.C., Quanser
Michel Lévis, M.A.S.C., Quanser

Quanser educational solutions
are powered by:



Course material
complies with:



CAPTIVATE. MOTIVATE. GRADUATE.

*ABET Inc. is the recognized accreditor for college and university programs in applied science, computing, engineering, and technology. ABET has provided leadership and quality assurance in higher education for over 75 years.

© 2012 Quanser Inc., All rights reserved.

Quanser Inc.
119 Spy Court
Markham, Ontario
L3R 5H6
Canada
info@quanser.com
Phone: 1-905-940-3575
Fax: 1-905-940-3576

Printed in Markham, Ontario.

For more information on the solutions Quanser Inc. offers, please visit the web site at:
<http://www.quanser.com>

This document and the software described in it are provided subject to a license agreement. Neither the software nor this document may be used or copied except as specified under the terms of that license agreement. All rights are reserved and no part may be reproduced, stored in a retrieval system or transmitted in any form or by any means, electronic, mechanical, photocopying, recording, or otherwise, without the prior written permission of Quanser Inc.

ACKNOWLEDGEMENTS

Quanser, Inc. would like to thank Dr. Hakan Gurocak, Washington State University Vancouver, USA, for his help to include embedded outcomes assessment.

CONTENTS

1	Introduction	4
2	Background	5
2.1	Modeling from First Principles	5
2.2	Desired Control Response	9
2.3	Ball and Beam Cascade Control Design	9
3	Pre-Lab Questions	14
4	Lab Experiments	15
4.1	Ideal PD Controller	15
4.2	Practical PD Controller	20
4.3	Results	24
5	System Requirements	25
5.1	Overview of Files	25
5.2	Software Setup	26
6	Lab Report	27
6.1	Template for Content (Ideal PD Control Experiments)	27
6.2	Template for Content (Practical PD Control Experiments)	28
6.3	Tips for Report Format	29

1 INTRODUCTION

The objective of this laboratory is to stabilize a ball to a desired position along a beam. Using the proportional-derivative (PD) family, a cascade control system is designed to meet a set of specifications.

Topics Covered

- Modeling dynamics of the ball from first-principles.
- Obtaining transfer function representation of the system.
- Design of a proportional-velocity (PV) compensator to control the position of the servo load shaft according to time-domain requirements.
- Assessment of how well the system meets design specifications using root locus.
- Design of a cascade control system to regulate the position of the ball on the beam.
- Simulation of the Ball and Beam control to ensure that the specifications are met without any actuator saturation.
- Implementation of the controllers on the Quanser BB01 device and evaluation of its performance.

Prerequisites

In order to successfully carry out this laboratory, the user should be familiar with the following:

- Data acquisition card (e.g., Q2-USB), the power amplifier (e.g., VoltPAQ-X1), and the main components of the SRV02 (e.g. actuator, sensors), as described in References [2], [3], and [1], respectively.
- Wiring and operating procedure of the SRV02 plant with the amplifier and DAQ devices, as discussed in Reference [1].
- Transfer function fundamentals, e.g. obtaining a transfer function from a differential equation.
- *SRV02 LabVIEW Integration* lab in [6] to get familiar with using **LabVIEW™** with the SRV02.

2 BACKGROUND

2.1 Modeling from First Principles

As illustrated in Figure 2.1, this system is comprised of two plants: the Rotary Servo (SRV02) and the Ball and Beam (BB01).




Figure 2.1: Ball and Beam open-loop block diagram.

The main objective in this section is to obtain the complete SRV02+BB01 transfer function

$$P(s) = P_{bb}(s) P_s(s) \quad (2.1)$$

where the BB01 transfer function is

$$P_{bb}(s) = \frac{X(s)}{\Theta_l(s)} \quad (2.2)$$

and the SRV02 transfer function is

$$P_s(s) = \frac{\Theta_l(s)}{V_m(s)}. \quad (2.3)$$

The BB01 transfer function describes the linear displacement of the ball, $x(t)$, with respect to the load angle of the servo, $\theta_l(t)$. In the next few sections, the time-based motion equations are developed and the transfer function is obtained. The SRV02 transfer function models the servo load gear position, $\theta_l(t)$, with respect to the servo input voltage, $V_m(t)$. Recall that in Modeling Laboratory ([6]), this transfer function was found to be:

$$P_s(s) = \frac{K}{s(\tau s + 1)} \quad (2.4)$$

The nominal model parameters, K and τ , of the SRV02 with no load and in high-gear configuration are:

$$K = 1.53 \text{ rad/(V.s)} \quad (2.5)$$

and

$$\tau = 0.0248 \text{ s} \quad (2.6)$$

Note: These parameters are different than the those found in the Modeling Laboratory because it does not include the inertial disc load.

2.1.1 Nonlinear Equation of Motion

In this section, the equation describing the motion of the ball, x , relative to the angle of the beam, α , is derived. Referring to Figure 2.2, the equation of motion, or *eom* for short, can be found starting from Newton's Law of Motion. The sum of forces acting on the ball along the beam equals:

$$m_b \left(\frac{d^2}{dt^2} x(t) \right) = \sum F \quad (2.7)$$




Figure 2.2: Free-body diagram of Ball and Beam.

where m_b is the mass of the ball.

Neglecting friction and viscous damping, the ball forces can be represented by

$$m_b \left(\frac{d^2}{dt^2} x(t) \right) = F_{x,t} - F_{x,r} \quad (2.8)$$

where $F_{x,r}$ is the force from the ball's inertia and $F_{x,t}$ is the translational force generated by gravity. For the ball to be stationary at a certain moment, i.e. be in equilibrium, the force from the ball's momentum must be equal to the force produced by gravity. As illustrated in Figure 2.3, the force $F_{x,t}$ in the x direction (along the beam) that is caused by gravity can be found as:




Figure 2.3: Completed Ball and Beam free-body diagram.

$$F_{x,t} = m_b g \sin \alpha(t) \quad (2.9)$$

The force caused by the rotation of the ball is

$$F_{x,r} = \frac{\tau_b}{r_b} \quad (2.10)$$

where r_b is the radius of the ball and τ_b is the torque which equals

$$\tau_b = J_b \left(\frac{d^2}{dt^2} \gamma_b(t) \right) \quad (2.11)$$

where γ_b is the ball angle. Using the sector formula, $x(t) = \gamma_b(t) r_b$, we can convert between linear and angular

displacement. Then, the force acting on the ball in the x direction from its momentum becomes:

$$F_{x,r} = \frac{J_b \left(\frac{d^2}{dt^2} x(t) \right)}{r_b^2} \quad (2.12)$$

Now, by substituting the rotational and translational forces into Equation 2.8, we can get the nonlinear equation of motion for the ball and beam as:

$$m_b \left(\frac{d^2}{dt^2} x(t) \right) = m_b g \sin \alpha(t) - \frac{J_b \left(\frac{d^2}{dt^2} x(t) \right)}{r_b^2} \quad (2.13)$$

Solving for the linear acceleration gives:

$$\frac{d^2}{dt^2} x(t) = \frac{m_b g \sin \alpha(t) r_b^2}{m_b r_b^2 + J_b} \quad (2.14)$$

2.1.2 Adding SRV02 Dynamics

In this section, the equation of motion representing the position of the ball relative to the angle of the SRV02 load gear is derived. The obtained equation will be nonlinear (includes a trigonometric term). Therefore, it will have to be linearized to use in control design.

Let's look at how we can find the relationship between the SRV02 load gear angle, θ_l , and the beam angle, α . Using the schematic given in *BB01 User Manual* ([4]), consider the beam and servo angles required to change the height of the beam by h . Taking the sine of the beam angle gives the expression

$$\sin \alpha(t) = \frac{h}{L_{beam}} \quad (2.15)$$

and taking the sine of the servo load shaft angle results in the equation

$$\sin \theta_l(t) = \frac{h}{r_{arm}} \quad (2.16)$$

From these we can obtain the following relationship between the beam and servo angle

$$\sin \alpha(t) = \frac{\sin \theta_l(t) r_{arm}}{L_{beam}} \quad (2.17)$$

To find the equation of motion that represent the ball's motion with respect to the SRV02 angle θ_l we need to linearize the equation of motion about the servo angle $\theta_l(t) = 0$. Insert the servo and beam angle relationship, Equation 2.17, into the nonlinear eom found in 2.14

$$\frac{d^2}{dt^2} x(t) = \frac{m_b g \sin \theta_l(t) r_{arm} r_b^2}{L_{beam} (m_b r_b^2 + J_b)} \quad (2.18)$$

About angle zero, the sine function can be approximated by

$$\sin \theta_l(t) \approx \theta_l(t) \quad (2.19)$$

Applying this to the nonlinear eom gives the *linear* equation of motion of the BB01

$$\frac{d^2}{dt^2} x(t) = \frac{m_b g \theta_l(t) r_{arm} r_b^2}{L_{beam} (m_b r_b^2 + J_b)} \quad (2.20)$$

To simplify the equation, we can lump the coefficient parameters of $\theta_l(t)$ into a single parameter K_{bb} . Here, K_{bb} is called the *model gain* of the Ball and Beam system.

2.1.3 Obtaining Transfer Function

In this section, we will derive the transfer function that describes the servo voltage to ball position displacement. The transfer function $P_{bb}(s)$ of the BB01 can be found by taking the Laplace transform of the linear equation of motion in 2.20 as:

$$P_{bb}(s) = \frac{X(s)}{\Theta_l(s)} = \frac{K_{bb}}{s^2} \quad (2.21)$$

As illustrated in Figure 2.1, both systems are in series. By inserting the BB01 plant transfer function 2.21 and the SRV02 plant transfer function 2.4 into Equation 2.1, we can derive the complete process transfer function $P(s)$ as:

$$P(s) = \frac{X(s)}{V_m(s)} = \frac{K_{bb} K}{s^3(\tau s + 1)} \quad (2.22)$$

This is the servo voltage to ball displacement transfer function.

Note: Finding the K_{bb} gain is left as an exercise in Section 3.

2.2 Desired Control Response

2.2.1 Time-Domain Specifications

The steady-state error, peak time, and percentage overshoot time-domain specifications for controlling the position of the SRV02 load shaft are:

$$e_{ss} = 0 \quad (2.23)$$

$$t_p = 0.15 \text{ s, and} \quad (2.24)$$

$$PO = 5.0 \% \quad (2.25)$$

Thus, when tracking the load shaft step reference, the transient response should have a peak time less than or equal to 0.15 seconds, an percentage overshoot (PO) less than or equal to 5 %, and no steady-state error.

The steady-state error, settling time, settling percentage, and percentage overshoot specifications for controlling the position of the ball are:

$$|e_{ss}| \leq 0.005 \text{ m} \quad (2.26)$$

$$t_s = 3.5 \text{ s} \quad (2.27)$$

$$c_{ts} = 0.04, \text{ and} \quad (2.28)$$

$$PO = 10.0 \% \quad (2.29)$$

Thus, given a step reference, the peak position of the ball (i.e., percentage overshoot) should not exceed 10%. After 3.5 seconds, the ball should settle within 4% of its steady-state value (i.e. not the reference) and the steady-state should be within 5 mm of the desired position.

2.3 Ball and Beam Cascade Control Design

The cascade control that will be used for the SRV02+BB01 system is illustrated by the block diagram in Figure 2.4. Based on the measured ball position $X(s)$, the ball and beam compensator, $C_{bb}(s)$ in the outer-loop computes the servo load shaft angle $\Theta_d(s)$ to attain the desired ball position $X_d(s)$. The inner loop is a servo position control system as described in the SRV02 Position Control Laboratory ([6]). Thus, the servo compensator $C_s(s)$ calculates the motor voltage required to track the desired load shaft angle.




Figure 2.4: Cascade control system used to control ball position in SRV02+BB01 plant.

2.3.1 Inner Loop Controller Design: SRV02 PV Position Controller

The inner loop implements a PV controller to manage the position of the SRV02 load shaft. In Figure 2.5 this is shown as $C_s(s)$. The design of this controller was explored previously in the SRV02 Position Control Laboratory. To meet the desired performance specifications given in Section 2.2.1, we need to compute the k_p and k_v (proportional and velocity) gains of the controller when the SRV02 is in the high-gear configuration.




Figure 2.5: SRV02 closed-loop system.

Recall that the standard definitions for peak time, t_p , and percent overshoot, PO , are:

$$t_p = \frac{\pi}{\omega_n \sqrt{1 - \zeta^2}} \quad (2.30)$$

and

$$PO = 100 \exp\left(\frac{-\pi \zeta}{\sqrt{1 - \zeta^2}}\right) \quad (2.31)$$

Using these formulas with the specifications in Equation 2.25, we can calculate the damping ratio needed to obtain the desired percent overshoot:

$$\zeta = 0.690 \quad (2.32)$$

And, the minimum natural frequency required to meet the desired peak time, given in 2.24, is

$$\omega_n = 28.9 \text{ rad/s} \quad (2.33)$$

We can find the necessary control gains by using the definitions for k_p and k_v from SRV02 Position Control Laboratory. The k_p gain is found from the model parameters in 2.5 and 2.6 as well as the desired natural frequency found in 2.33:

$$k_p = 13.5 \text{ V/rad} \quad (2.34)$$

The velocity gain k_v is obtained using the model parameters and the minimum damping ratio specification in 2.32:

$$k_v = 0.078 \text{ V.s/rad} \quad (2.35)$$

2.3.2 Outer Loop Controller Design

We plan to use a dynamic compensator in the outer loop in Figure 2.6. Generally speaking, adding a zero in the forward-path increases the bandwidth of the closed-loop system. Adding a pole increases the rise time and overshoot of the system and makes it overall less stable. In our case, the bandwidth must be increased and the overshoot has to be minimized.

We can simplify the cascade controller shown in Figure 2.4, by replacing the inner loop with a single block $G_s(s)$ as:

For now, we will consider the dynamics of the inner loop as negligible ($G_s(s) = 1$). Therefore, it is assumed that the inner loop can control the position of the SRV02 load shaft perfectly, hence the desired load shaft angle equals the actual load shaft angle:

$$\theta_l(t) = \theta_d(t) \quad (2.36)$$

We will design a PD controller for the "BB01 Compensator" block in Figure 2.6. This controller is in the outer loop. Given the desired ball position, $X_d(s)$, and the actual ball position $X(s)$, the error in the system will be used by the controller to generate $\Theta_d(s)$ as a desired input for the inner loop controller $G_s(s)$. Since we assumed that the inner loop controller can perfectly track the desired shaft angle, it will manage the dynamics of the SRV02 and move the SRV02 shaft to the $\Theta_l(s)$ position.




Figure 2.6: BB01 closed-loop system.

A *traditional* PD controller is given by $G_{PD}(s) = K_c \cdot (s + z)$. However, before we can proceed, we need to make a slight variation in the traditional PD controller. As we have done previously in the SRV02 Position Control Laboratory, we will replace the traditional PD controller with the one shown in Figure 2.7. Since this version does not directly feed the setpoint velocity $s X_d(s)$, it is more desirable in practical applications due to the noise in the signals.




Figure 2.7: Ball and beam ideal PD controller (Also called PV controller).

Ideal PD Control Design

Using the block diagram, the outer-loop controller can be written as:

$$\Theta_d(s) = K_c (z (X_d(s) - X(s)) - s X(s)) \quad (2.37)$$

Recall that $\Theta_l(t) = \Theta_d(t)$ when the servo dynamics is ignored. By substituting the above transfer function into 2.21 and solving for $X(s)/X_d(s)$, we can obtain the BB01 closed-loop transfer function as:

$$\frac{X(s)}{X_d(s)} = \frac{K_{bb} K_c z}{s^2 + K_{bb} K_c s + K_{bb} K_c z} \quad (2.38)$$

In Figure 2.7, z and K_c are the zero location and the gain for the controller, respectively. We can compute the value for z to meet the settling time and overshoot specifications. Also, we can find the value of K_c to satisfy a given natural frequency and damping ratio.

A standard second order system is given by the following transfer function:

$$\frac{Y(s)}{R(s)} = \frac{\omega_n^2}{s^2 + 2\zeta\omega_n s + \omega_n^2} \quad (2.39)$$

By matching the coefficients of the denominator of this equation to those in equation 2.38, we can find the expressions for z and K_c to meet a certain ω_n and ζ :

$$z = \frac{\omega_n^2}{K_{bb} K_c} \quad (2.40)$$

and the gain must be

$$K_c = \frac{2\zeta\omega_n}{K_{bb}} \quad (2.41)$$

Practical PD Controller

and

$$p_c = -\omega_f \quad (2.49)$$

respectively. Since $p_c < z_c$, the pole comes before the zero along the negative real-axis, this is a *lead* type compensator.

A standard third-order characteristic equation can be written as:

$$(s^2 + 2\zeta\omega_n s + \omega_n^2)(1 + T_p s) \quad (2.50)$$

where T_p is the pole decay in seconds. When expanded, the third-order characteristic equation becomes

$$s^3 + \frac{(2\zeta\omega_n T_p + 1)}{T_p} s^2 + \frac{(\omega_n^2 T_p + 2\zeta\omega_n)}{T_p} s + \frac{\omega_n^2}{T_p} \quad (2.51)$$

The characteristic equation of the closed-loop system (equation 2.45) using the *practical* PD controller is

$$s^3 + s^2\omega_f + (K_{bb}K_c\omega_f + K_{bb}K_c z) s + K_{bb}K_c z\omega_f \quad (2.52)$$

The following expression is obtained after equating the coefficients of the s^2 terms in equations 2.51 and 2.52,

$$\omega_f = \frac{2\zeta\omega_n T_p + 1}{T_p} \quad (2.53)$$

Solving for T_p gives the the required pole location given the desired response specifications

$$T_p = \frac{1}{\omega_f - 2\zeta\omega_n} \quad (2.54)$$

3 PRE-LAB QUESTIONS

- Find K_{bb} by simplifying the expression given in Equation 2.18. Then, evaluate it using the system parameters given in [4]. **Hint:** Recall that the mass moment of inertia of a solid sphere is $J = \frac{2mr^2}{5}$.
- Find the steady-state error of the the Ball and Beam system given by the $P_{bb}(s)$ transfer function. The system is shown in Figure 3.1. The compensator is unity




Figure 3.1: Unity feedback system.

$$C(s) = 1 \quad (3.1)$$

and the reference step is

$$R(s) = \frac{R_0}{s} \quad (3.2)$$

where R_0 is the step amplitude. Note that in this calculation the SRV02 dynamics is to be ignored and only the BB01 plant is to be considered.

- Using Figure 2.6, find the closed-loop transfer function of the BB01 system with proportional control $C_{bb}(s) = K_c$.
- Plot the root locus of the BB01 plant $P_{bb}(s)$. Describe how the poles behave as K_c goes to infinity.
- Find the natural frequency and damping ratio required to achieve the time-domain specifications of the Ball and Beam plant given in Section 2.2.1.
- After you plot the root locus of the BB01 plant $P_{bb}(s)$, describe where the poles should be to satisfy the desired response specifications.
- Discuss the response if the poles lie beyond the radius circle along the diagonal lines, i.e. away from the imaginary axis. Also, comment on what happens if the poles of the system lie inside the diagonal lines along the radius circle, i.e. moving towards the real axis. Make references to its effects on the settling time and overshoot of the response.
- Based on the root locus obtained in question 6 previously, can the specifications of the Ball and Beam system be satisfied using a proportional controller? Discuss.
- Assume that a *traditional* PD controller ($C_{bb}(s) = K_c(s+z)$) is used in the system given in Figure 2.6. Find the BB01 error transfer function.
- Find the steady-state error of the BB01 closed-loop system with the *traditional* PD controller. Can the steady-state error requirement in 2.26 be satisfied?
- Based on the expressions found in Equations 2.40 and 2.41, evaluate numerically the zero location and gain needed to satisfy the specifications.
- Find expressions for the zero location, z , and the compensator gain, K_c for the *practical* PD controller to satisfy ω_n and ζ in Section 2.2.1 and the desired filter cutoff frequency in 2.43. Then, evaluate numerically the pole time constant, zero location and gain needed to satisfy the specifications.

4 LAB EXPERIMENTS

The main goal of this laboratory is to explore position control of the Ball and Beam system using cascade control.

In this laboratory, you will conduct two groups of experiments:

1. *Ideal* PD controller
 - simulation with no servo dynamics,
 - simulation with servo dynamics
2. *Practical* PD controller and servo dynamics
 - simulation with servo dynamics,
 - implementation

4.1 Ideal PD Controller

4.1.1 Simulation with No Servo Dynamics

Root Locus of Ideal Control

The *BB01 Root Locus - Ideal (Student)* VI, shown in Figure 4.1 plots the root locus of the forward path when using the *ideal PD* control.




Figure 4.1: VI used to plot root locus of BB01 system with ideal control.

1. In the *Ball and Beam Control (Student)* LabVIEW project, open the *BB01 Root Locus - Ideal (Student)* VI.
2. Enter the SRV02 model parameters K and τ found in the SRV02 Modeling lab (see [6]), or just use the loaded nominal parameters.
3. Enter the BB01 model gain found in Pre-Lab question 1 in the K_{bb} control box on the front panel.
4. Enter the BB01 compensator gain K_c and the compensator zero, z , that were found in Pre-Lab question 11.

- Using **LabVIEW™**, plot the root locus of BB01 loop transfer function $L(s) = C_{bb}(s)P_{bb}(s)$ when using the *ideal PD* controller. As illustrated in Figure 4.2 the MathScript node is not complete, i.e., the plant and controller transfer functions are integrators. Enter the BB01 plant in the P_{bb} variable and the ideal PD compensator in the C_{bb} variable. In this case, use the ideal compensator $C_{bb}(s) = K_c(s + z)$. In the root locus plot, ensure the poles go through the desired locations at the gain that was computed in Pre-Lab question 11.




Figure 4.2: MathScript used to generate ideal loop transfer function of BB01 system (incomplete).

Experimental Setup for Simulation

The *BB01 Ideal Control Simulation VI* shown in Figure 4.3 is used to simulate the closed-loop position response of the Ball and Beam when using the outer-loop control. In this case the SRV02 dynamics are neglected, i.e. $\theta_d = \theta_l$. The response is simulated using the developed nonlinear model of the Ball and Beam.




Figure 4.3: VI used to simulate the outer closed-loop BB01 system.

The *BB01 Model* sub-VI includes the $P_{bb}(s)$ transfer function that was derived in Section 2.1.3. Recall that in Section 2.1.2 the model had to be linearized in order to obtain the $P_{bb}(s)$ transfer function. This nonlinearity is re-introduced in the *BB01 Model* subsystem in order to represent the plant more accurately and ensure the specifications can still be satisfied. The *ideal PD* compensator designed in Section 2.3.2 is also implemented. Note that it includes a *Saturation* block that limits the SRV02 angle between ± 56 degrees.

IMPORTANT: Before you can conduct these experiments, you need to make sure that the lab files are configured according to your system setup. If they have not been configured already, then you need to go to Section 5.2 to configure the lab files first.

Simulation

The closed-loop step position response of the BB01 is simulated to verify that the specifications are met. As previously mentioned, the simulation is performed using the nonlinear model of the Ball and Beam and the *ideal PD* controller that was designed.

1. In the *Ball and Beam Control (Student)* LabVIEW project, open the *BB01 Ideal Control Simulation VI* shown in Figure 4.3. Make sure it is setup as described in Section 5.2.
2. Enter the SRV02 model parameters K and τ found in the SRV02 Modeling lab (see [6]), or just use the loaded nominal parameters.
3. Enter the BB01 model gain found in Pre-Lab question 1 in the K_{bb} control box on the front panel.
4. Enter the BB01 compensator gain K_c and the compensator zero, z , that were found in Pre-Lab question 11.
5. To set the desired ball position to a ± 5 cm square wave at 0.05 Hz, set the *Signal Generator* controls to the following:
 - *Signal Type* = square
 - *Amplitude* = 5 cm
 - *Offset* = 0
 - *Frequency* = 0.05 Hz
6. Run the VI. Using the default control parameters, the scopes should be displaying responses similar to Figure 4.3. In the *Ball Position* scope, the blue trace is the desired ball position and the red trace is the simulated response. The bottom chart plots the SRV02 angular position.
7. Click on the STOP button when you have obtained a suitable response.
8. Generate figures showing the nonlinear Outer-Loop BB01 ball position response and the corresponding servo angle. To save the response of a plot after the VI is stopped, right-click on a chart, go to Export | Export Simplified Image, and save the image as a BMP on the clipboard. This can then be saved in your graphics editor or to your report.
9. Measure the steady-state error, the settling time, and the percent overshoot of the simulated response. Use the cursors in the *Ball Position (cm)* graph in the *Measure* tab to take your measurements.
10. Does the outer-loop ideal PD response satisfy the specifications given in Section 2.2.1 while keeping the servo angle between ± 56 degrees? Some tolerance is allowed on the settling time specification: it should not exceed 3.75 seconds (rather than 3.50 seconds). If the steady-state error and percent overshoot do not meet the desired specifications and the settling time goes over 3.75 seconds then go back to your control design. If the settling time does not satisfy the original specifications but is kept below the allowed tolerance, explain any possible source for this discrepancy.

4.1.2 Simulation with Servo Dynamics

Experimental Setup

The servo dynamics can now be added and the closed-loop position response with the cascade control system can be simulated using the LabVIEW *SRV02+BB01 Control Simulation VI* pictured in Figure 4.4. This VI is for the block diagram shown in Figure 2.4.




Figure 4.4: VI used to simulate cascade control system.

This VI includes the nonlinear model of the BB01 plant and the transfer function representing the SRV02 voltage-to-position relationship. The proportional-velocity position controller designed in Section 2.3.1 is implemented in the *SRV02 PV Position Control* block.

The cascade controller is the algorithm that will be implemented on the actual SRV02+BB01 device. Before deployment, we need to confirm that the specifications are still satisfied when the servo dynamics are added. In addition, the servo angle must be kept between ± 56 degrees and the servo voltage cannot exceed ± 10 V.

Simulation

The purpose of this simulation is to see how the settling time, overshoot, and steady-state error of the response changes when the inner-loop servo control that was designed in Section 2.3.1 is added. In this simulation, you will be using the *ideal PD* controller.

1. In the *Ball and Beam Control (Student)* LabVIEW project, open the *SRV02+BB01 Control Simulation* VI. Make sure it is setup as described in Section 5.2.
2. Enter the SRV02 model parameters K and τ found in the SRV02 Modeling lab (see [6]), or just use the loaded nominal parameters.
3. Enter the BB01 model gain found in Pre-Lab question 1 in the K_{bb} control box on the front panel of the VI.
4. Enter the BB01 compensator gain K_c and the compensator zero, z , that were found in Pre-Lab question 11.
5. Enter the SRV02 PV gains: called variables K_p and K_v in from Section 2.3.1.
6. To set the desired ball position to a ± 5 cm square wave at 0.05 Hz, set the *Signal Generator* controls to the following:
 - *Signal Type* = square
 - *Amplitude* = 5 cm
 - *Offset* = 0
 - *Frequency* = 0.05 Hz

7. Set the *Select Control* switch to the downward position in order to use the *ideal PD* controller when simulating.
8. Run the VI. The scopes should be displaying responses similar to Figure 4.4. In the *Ball Position* scope, the blue trace is the desired ball position and the red trace is the simulated response. The bottom left chart plots the desired (blue) and simulated (red) SRV02 angular positions and the bottom right chart displays the DC motor voltage.
9. Click on the STOP button when you have obtained a suitable response.
10. Plot the *Ideal PD* cascade ball position, servo angle, and servo input voltage response.
11. Measure the steady-state error, the settling time, and the percent overshoot of the ideal PD cascade control response. Use the cursors in the *Ball Position (cm)* graph in the *Measure* tab to take your measurements.
12. Are the specifications in Section 2.2.1 still satisfied after adding the servo dynamics? Also, make sure the servo angle is within ± 56.0 degrees and the servo voltage is between ± 10.0 V. *Do not go back in the control design if some specifications are not met.*

4.2 Practical PD Controller

4.2.1 Simulation with Servo Dynamics

Root Locus

Follow these steps to generate the root locus of the practical control.

1. In the *Ball and Beam Control (Student)* LabVIEW project, open the *BB01 Root Locus - Practical (Student)* VI. Make sure it is setup as described in Section 5.2.
2. As detailed in Section 4.1.2, make sure the SRV02 modeling parameters K and τ , the BB01 model gain K_{bb} , and the BB01 control gains K_c and z are entered in the VI.
3. The filter cutoff, ω_f , should already be set as $\omega_f = 2\pi = 6.28$ rad/s in the VI (See Section 5.2 for more details).
4. Use the VI to plot the root locus of BB01 loop transfer function when using the *practical PD* compensator. As illustrated in Figure 4.2 the MathScript node is not complete, i.e., the plant and controller transfer functions are set to integrators. Enter the BB01 plant in the P_{bb} variable and the *practical* controller in the C_{bb} variable. The practical PD compensator, $C_{bb}(s)$, is given in Equation 2.47. Show the desired locations of the poles on the plot and ensure the poles go through the desired locations at the gain that was computed.

Simulation

The *practical* PD controller developed in Section 2.3.2 is simulated in this section. This is the compensator that will be used to control the actual BB01 device. The control gain and zero may have to be fine-tuned in order to compensate for the added dynamics of the filtering and the inner-loop servo control. Follow these steps to simulate the closed-loop *practical cascade PD* response:

1. In the *Ball and Beam Control (Student)* LabVIEW project, open the *SRV02+BB01 Control Simulation* VI. Make sure it is setup as described in Section 5.2.
2. As in detailed in Section 4.1.2, make sure the SRV02 modeling parameters K and τ , the BB01 model gain K_{bb} , the BB01 control gains K_c and z , and the SRV02 PV gains K_p and K_v are entered in the VI.
3. The filter cutoff, ω_f , should already be set as $\omega_f = 2\pi = 6.28$ rad/s in the VI (See Section 5.2 for more details).
4. Set the desired ball position to a ± 5 cm square wave at 0.05 Hz (as in Section 4.1.2).
5. To simulate using the *practical PD* controller, set the *Select Control* switch to the upward position.
6. Run the VI. The VI is shown running in Figure 4.4 (when using the ideal controller).
7. Plot the *practical cascade* ball position, servo angle, and servo input voltage response (use the *Export* command).
8. Measure the steady-state error, the settling time, and the percent overshoot of the simulated *practical cascade PD control* response.
9. Does the simulated response satisfy the specifications given in Section 2.2.1 while keeping the servo angle between ± 56.0 degrees and the servo voltage between ± 10.0 V?
10. If the specifications were satisfied, then jump to Step 12 below. If the requirements are not met, then the controller needs to be tuned. One method is to redesign the compensator gain, zero location, and pole time constant, according to more stringent restrictions. For instance, try simulating the system for a K_c , z , and T_p , generated according to a percent overshoot of 8% instead of 10%. To do this, open the *BB01 Calc Practical Control Gains (Student)* VI and complete the MathScript to compute the gains automatically according to a given set of percent overshoot, settling time, and filter cutoff frequency specifications. The incomplete script

is shown below. Enter the appropriate equations for T_p , z , and K_c . Then, simulate the system and see if the specifications are satisfied. Note that the cutoff filter frequency should remain as specified in Equation 2.43, or 1 Hz.

```
% BB01 high-pass filter cut-off frequency (rad/s)
wf = 2*pi*1;
% Damping ratio from overshoot specification
zeta = -log(P0/100) * sqrt( 1 / ( ( log(P0/100) )^2 + pi^2 ) );
% Natural frequency from specifications (rad/s)}
wn = -log( c_ts * (1-zeta^2)^(1/2) ) / (zeta * ts);
% Resulting pole decay (s)
Tp = 1;
% Zero location (rad/s)
z = 1;
% Compensator Gain (rad/m)
Kc = 1;
```

11. **Tuned Practical PD (Simulation):** Record the gain and zero that have been fine-tuned for the response to meet the specifications along with the new specifications used to generate those control parameters.
12. Plot the simulated response.
13. Give the resulting steady-state error, settling time, and percent overshoot of the response.

4.2.2 Implementation with Practical PD Controller

In this experiment, the position of the ball on the BB01 device will be controlled using the developed *practical* PD controller. Measurements will then be taken to ensure that the specifications are satisfied.

Experimental Setup

The *BB01 Control VI* shown in Figure 4.5 is used to perform the position control exercises in this laboratory. It includes drivers that interface with the DC motor and sensors of the Ball and Beam system. The VI also implements *practical PD* control detailed in Section 2.3.2.

IMPORTANT: Before you can conduct this experiment, you need to make sure that the lab files are configured according to your system setup. If they have not been configured already, then you need to go to Section 5.2 to configure the lab files first.

1. In the *Ball and Beam Control (Student)* LabVIEW project, open the *SRV02+BB01 Control Simulation VI*. Make sure it is setup as described in Section 5.2.
2. As detailed in Section 4.1.2, make sure the SRV02 modeling parameters K and τ , the BB01 model gain K_{bb} , the BB01 control gains K_c and z , and the SRV02 PV gains K_p and K_v are entered in the VI.
3. The filter cutoff, ω_f , should already be set as $wf = 2\pi = 6.28$ rad/s in the VI (See Section 5.2 for more details).
4. Set the desired ball position to a ± 5 cm square wave at 0.05 Hz (as in in Section 4.1.2).
5. Place the *Ball Command* switch to the upward position in order to generate the setpoint using the VI signal generator.
6. Turn ON the power amplifier.
7. Run the VI. The scopes should be displaying responses similar to Figure 4.5.
8. When a suitable response is obtained, click on the *STOP* button to stop running the code. Plot the ball position and servo angle response as well as the input voltage.




Figure 4.5: VI that runs the practical PD controller on the Ball and Beam system.

9. **Tuned Practical PD (Implementation):** Measure the steady-state error, the settling time, and the percent overshoot. Does the response satisfy the specifications given in Section 2.2.1? Give one reason why the designed gain and zero could fail to give a successful closed-loop response on the actual system?
10. If the specifications have been satisfied, then go to the next step. Otherwise tune the gains as described in Section 4.2.1 and run the experiment again until a satisfactory response is obtained. In the case where the steady-state error is not satisfied, integral action can be introduced in the outer-loop controller. To do this, increase the ki control box on the VI front panel in small increments until the error is minimized enough. Briefly explain the procedure to get those new control parameters (including the integral gain, if necessary) and give the gain and zero used to obtain the response.
11. Make sure the VI is stopped.
12. Shut off the power amplifier if no more experiments will be performed in this session.

4.2.3 Controller using the Remote Sensor (Optional)

In this experiment, the position of the ball on the BB01 device will be controlled using the developed *practical PD control* but the setpoint is given with the remote sensor. This experiment can be conducted only if the remote sensor SS01 module detailed in [4] is available.

1. Follow steps 1-4 in Section 4.2.2 to setup the *BB01 Control VI*.
2. Place the *Ball Command* switch to the DOWN position in order to generate the setpoint using the remote sensor (SS01) module.
3. Run the VI.
4. Move the ball back and forth on the remote sensor and observe the response obtained in the scopes. Figure 4.6 shows the VI running and a sample response.
5. When done, click on the *STOP* button.
6. Shut off the power of the amplifier if no more experiments will be performed in this session.




Figure 4.6: BB01 ball position response with SS01.

4.3 Results

Fill out Table 4.1 below with your answers to the Pre-Lab questions and your results from the lab experiments.

Section / Question	Description	Symbol	Value	Unit
Question 11	Pre-Lab: Ideal PD Control Design Compensator Gain Compensator Zero	K_c Z		
Question 12	Pre-Lab: Practical PD Control Design Compensator Gain Compensator Zero Compensator Pole Time Constant	K_c Z T_p		
4.1.1	In-Lab Simulation: Cascade ideal PD with no servo dyn. Steady-state error Settling time Percentage overshoot	e_{ss} t_s PO		
4.1.2	In-Lab Simulation: Cascade ideal PD with servo dyn. Steady-state error Settling time Percentage overshoot	e_{ss} t_s PO		
4.2.1	In-Lab Simulation: Cascade Practical PD Steady-state error Settling time Percentage overshoot	e_{ss} t_s PO		
Section 4.2.1, step 11	In-Lab Simulation: Cascade Tuned Practical PD Compensator Gain Compensator Zero Steady-state error Settling time Percentage overshoot	K_c Z e_{ss} t_s PO		
Section 4.2.2, step 9	In-Lab Implementation: Tuned Practical PD Steady-state error Settling time Percentage overshoot	e_{ss} t_s PO		

Table 4.1: Summary of results for the Ball and Beam position laboratory.

5 SYSTEM REQUIREMENTS

Required Hardware

- Data-acquisition (DAQ) device that is compatible with LabVIEW™ and Quanser Rapid Control Prototyping Toolkit®.
- Quanser SRV02-ET rotary servo. See Reference [1].
- Quanser VoltPAQ power amplifier, or equivalent (e.g. Reference [3] for VoltPAQ User Manual).

Required Software

- NI LabVIEW™
- NI LabVIEW Control Design and Simulation Module
- NI LabVIEW MathScript RT Module or later
- For NI CompactRIO users:
 - NI LabVIEW Real-Time Module
 - NI LabVIEW FPGA Module
 - NI-RIO Drivers
- Quanser Rapid Control Prototyping Toolkit®

5.1 Overview of Files

File Name	Description
Ball and Beam User Manual.pdf	This manual describes the hardware of the Ball and Beam system and explains how to setup and wire the system for the experiments.
Ball and Beam Control (STUDENT).pdf	This laboratory guide contains pre-lab questions and lab experiments demonstrating how to design and implement a position controller on the Quanser SRV02 Ball and Beam plant using LabVIEW™.
BB01 Root Locus - Ideal (Student).vi	Plots the root locus of the ideal PD control forward path. Must be completed by student.

(Continued on the next page)

File Name	Description
BB01 Ideal Control Simulation.vi	Simulates the closed-loop system when using only the outer-loop ball position controller with the BB01 system, i.e. no inner loop control of the servo position.
BB01 Root Locus - Practical (Student).vi	Plots the root locus of the practical PD control forward path. Must be completed by student.
BB01 Calc Practical Controls Gain (Student).vi	Computes the practical PD control parameters based on the specifications give. Must be completed by student.
SRV02+BB01 Control Simulation.vi	Simulates the cascade ball position controller. Both the outer-loop ball position control and the inner-loop servo position control are used in this file.
BB01 Control	Implements the closed-loop cascade position controller on the actual BB01 system using LabVIEW™.

Table 5.1: Files supplied with the SRV02 Speed Control laboratory.

5.2 Software Setup

Follow these steps to get the system ready for this lab:

1. Load the LabVIEW™ software.
2. Open the LabVIEW project called *Ball and Beam Control (Student).lvproj* shown in Figure 5.1.
3. To run the Ball and Beam experiments, open the *BB01 Control VI* under My Computer.
4. **Configure DAQ Device:** Before running the VI, make sure you set the correct *Board type* (e.g., 'q1_cRIO', 'q2_usb', 'q8_usb', 'qpid', or 'qpid_e') in the HIL Initialize block.
5. **NI CompactRIO Users:** Before running the VI, make sure you can connect to your CompactRIO through the Measurement & Automation software. See the SRV02 cRIO User Manual ([5]).
6. **Channel Configuration:** For any of these VIs, the analog input and output channels are set, by default, to match the wiring in the Ball and Beam User Manual ([4]). If the wiring is different on your system, make sure the VI uses the correct channels. For instance, if your ball sensor is connected to Analog Input Channel #1 on your DAQ, then set the *Ball Position Channel* control box in the VI to 1 (instead of 0).




Figure 5.1: Ball and Beam Project.

6 LAB REPORT

This laboratory contains two groups of experiments, namely,

1. *Ideal* PD controller with and without servo dynamics (outer-loop only and cascade control), and
2. *Practical* PD controller with servo dynamics (cascade control).

For each experiment, follow the outline corresponding to that experiment to build the *content* of your report. Also, in Section 6.3 you can find some basic tips for the *format* of your report.

6.1 Template for Content (Ideal PD Control Experiments)

I. PROCEDURE

1. *Simulation with no servo dynamics*
 - Briefly describe the main goal of the simulation.
 - Briefly describe the simulation procedure (Section 4.1.1)
2. *Simulation with servo dynamics*
 - Briefly describe the main goal of this simulation.
 - Briefly describe the simulation procedure (Section 4.1.2)

II. RESULTS

Do not interpret or analyze the data in this section. Just provide the results.

1. Response plot from step 8 in Section 4.1.1, *Simulation with no servo dynamics*
2. Response plot from step 10 in Section 4.1.2, *Simulation with servo dynamics*
3. Provide applicable data collected in this laboratory (from Table 4.1).

III. ANALYSIS

Provide details of your calculations (methods used) for analysis for each of the following:

1. Root locus plot from step 5 in Section 4.1.1, *Simulation with no servo dynamics*
2. Steady-state error, settling time and percent overshoot in step 9 in Section 4.1.1, *Simulation with no servo dynamics*
3. Steady-state error, settling time and percent overshoot in step 11 in Section 4.1.2, *Simulation with servo dynamics*

IV. CONCLUSIONS

Interpret your results to arrive at logical conclusions for the following:

1. Whether the controller meets the specifications in step 10 in Section 4.1.1, *Simulation with no servo dynamics*
2. Whether the controller meets the specifications in step 12 in Section 4.1.2, *Simulation with servo dynamics*

6.2 Template for Content (Practical PD Control Experiments)

I. PROCEDURE

1. *Simulation*

- Briefly describe the main goal of the simulation.
- Briefly describe the simulation procedure (Section 4.2.1)

2. *Implementation*

- Briefly describe the main goal of this experiment.
- Briefly describe the experimental procedure (Section 4.2.2)

II. RESULTS

Do not interpret or analyze the data in this section. Just provide the results.

1. Response plot from step 7 in Section 4.2.1, *Simulation with practical PD controller*
2. Response plot from step 12 in Section 4.2.1, for *Tuned Practical PD (Simulation)*
3. Response plot from step 8 in Section 4.2.2, *Implementation with tuned practical PD controller*
4. Provide applicable data collected in this laboratory (from Table 4.1).

III. ANALYSIS

Provide details of your calculations (methods used) for analysis for each of the following:

1. Steady state error, the settling time and percent overshoot in step 8 in Section 4.2.1, *Simulation with practical PD controller*
2. Steady state error, the settling time and percent overshoot in step 13 in Section 4.2.1, *Simulation with practical PD controller*
3. Steady state error, the settling time and percent overshoot in step 9 in Section 4.2.2, *Implementation with practical PD controller*

IV. CONCLUSIONS

Interpret your results to arrive at logical conclusions for the following:

1. Whether the controller meets the specifications in step 9 in Section 4.2.1, *Simulation with practical PD controller*
2. Whether the controller meets the specifications in step 9 in Section 4.2.2, *Implementation with practical PD controller*

6.3 Tips for Report Format

PROFESSIONAL APPEARANCE

- Has cover page with all necessary details (title, course, student name(s), etc.)
- Each of the required sections is completed (Procedure, Results, Analysis and Conclusions).
- Typed.
- All grammar/spelling correct.
- Report layout is neat.
- Does not exceed specified maximum page limit, if any.
- Pages are numbered.
- Equations are consecutively numbered.
- Figures are numbered, axes have labels, each figure has a descriptive caption.
- Tables are numbered, they include labels, each table has a descriptive caption.
- Data are presented in a useful format (graphs, numerical, table, charts, diagrams).
- No hand drawn sketches/diagrams.
- References are cited using correct format.

REFERENCES

- [1] Quanser Inc. Srv02 user manual. 2009.
- [2] Quanser Inc. Q2-usb User's Manual. 2010.
- [3] Quanser Inc. Voltpaq User's Manual. 2010.
- [4] Quanser Inc. Ball and beam user manual. 2011.
- [5] Quanser Inc. *SRV02 cRIO User Manual*, 2011.
- [6] Quanser Inc. SRV02 lab manual. 2011.

Over ten rotary experiments for teaching fundamental and advanced controls concepts



Quanser's rotary collection allows you to create experiments of varying complexity – from basic to advanced. Your lab starts with the Rotary Servo Base Unit and is designed to help engineering educators reach a new level of efficiency and effectiveness in teaching controls in virtually every engineering discipline including electrical, computer, mechanical, aerospace, civil, robotics and mechatronics. For more information please contact info@quanser.com

©2013 Quanser Inc. All rights reserved.



INFO@QUANSER.COM

+1-905-940-3575

QUANSER.COM

Solutions for teaching and research. Made in Canada.